# Canada Communicable Disease Report

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Vol. 25-24

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## Contained in this issue:

Measles Surveillance: Guidelines for Laboratory Support . . . . . 201

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### MEASLES SURVEILLANCE: GUIDELINES FOR LABORATORY SUPPORT

#### Introduction

Following the commitment by all provinces and territories to eliminate measles by the year 2005, a two-dose measles vaccine schedule has been adopted across the country, augmented by mass vaccination "catch-up" programs in almost all jurisdictions. As the number of cases of measles declines, the importance of surveillance will become even greater than at present. It will be more and more crucial that all suspected cases of measles be reported and samples from all sporadic cases be submitted for full laboratory investigation.

To oversee and document the process of measles elimination, the Laboratory Centre for Disease Control (LCDC) convened the Working Group on Measles Elimination in Canada (WGMEC) in 1995. This group reviews measles elimination efforts across the country, monitors case incidence, and makes recommendations about national surveillance.

The importance of a standard set of laboratory procedures to be used across the country is evident in any assessment of measles surveillance procedures. A subgroup of WGMEC\*, with additional laboratory support, has therefore developed these guidelines to ensure optimum laboratory surveillance of measles at a time when, with fewer cases, false positive serological results will be more likely, and it will become increasingly important to detect any importation of virus into a community.

Effective laboratory support for surveillance requires that

- health-care professionals be aware of the system and its requirements,
- appropriate specimens are collected and sent to a laboratory capable of performing the necessary tests,
- reliable equipment is used and tests performed correctly,

\* The members of the Working Group on Measles Elimination in Canada laboratory issues subgroup are: Dr. M. Fearon (Ontario Ministry of Health), Dr. B. Ward (McGill University), Dr. S. Ratnam (Newfoundland Public Health Laboratory), Dr. G. Tipples (LCDC), M. Fauvel (Laboratoire de santé publique du Québec), Dr. J. Waters (Alberta Health), and M. Litt (LCDC).

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- feedback is provided to the appropriate authorities and in a timely fashion, and
- the integrity of the surveillance system is monitored on an ongoing basis.

Both the public-health authority and the laboratory have a role in ensuring optimum laboratory surveillance of measles in a jurisdiction. These guidelines have been structured to reflect these roles. The role of public-health professionals is to educate health-care providers on the importance of laboratory testing and to make sure that the laboratory resources are accessible. The role of laboratories is to develop and maintain testing resources, standard procedures, and quality assurance; and to analyze data and provide feedback of the results to public-health authorities.

**The following guidelines, which have been updated in a number of key areas to reflect the progress of measles elimination in Canada, propose a standard set of responsibilities for adoption in each province or territory to facilitate effective laboratory support for the surveillance of measles.**

## Laboratory Issues

### *Laboratory Methods for Diagnosing Measles*

- Laboratory diagnosis of measles is done by serologic testing and/or culture.
- Measles specific IgM serology is the standard test of choice for routine diagnosis of measles.
- Demonstration of a significant increase\* in measles specific IgG titre is a reliable alternative serologic method for diagnosis.

Measles virus isolates are important for surveillance purposes (molecular epidemiology) as well as for confirmatory purposes; therefore, collection of appropriate specimens for culture is recommended for all sporadic cases and a sample of cases in an outbreak.

Appendix A contains the measles case definitions and appendix B outlines the decision-making process for measles case investigation.

All clinical and suspected cases of measles should be confirmed by a laboratory. If an epidemiologic link to an already confirmed case has been established, laboratory confirmation is not necessary for the case to meet the “confirmed” case definition. To detect cases of measles and to avoid misdiagnosis of rash illness, it is important to use a laboratory testing protocol to screen for common red rash illnesses.

**Upon suspicion of measles, clinicians should immediately collect specimens for serology and virus isolation for the purpose of laboratory confirmation.**

\* A four-fold rise in IgG titre using such assays as complement fixation (CF) or hemagglutination (HI) has traditionally been used to define significant increases in antibody titres. These tests have mostly been replaced by enzyme immunoassays (EIAs), which may or may not permit quantitative measurement. Many commercial IgG EIAs, such as the Behring Enzygnost test, use a slightly modified procedure to demonstrate a significant rise in antibody titre between acute and convalescent serum specimens. Alternatively, paired specimens can be referred to a reference center for CF, HI, or plaque reduction neutralization testing.

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### *Information on Measles Serologic Testing*

- Measles specific IgM antibodies appear around the time of rash onset and persist for at least 28 days.
- Measles specific IgG antibodies appear around the same time as IgM antibodies.
- A single blood specimen collected 3 to 28 days after rash onset is usually satisfactory for IgM serology.
- IgG serology using paired “acute” and “convalescent” specimens is a reliable test for diagnosing measles, provided that specimens are collected at the appropriate times.
- For IgG serology, the first (acute) sample should be obtained as soon as possible after the onset of the rash and, in any event, no later than 7 days afterwards. The second (convalescent) sample should be collected 10 to 20 days after the first sample. These paired sera must be tested simultaneously.

The timing of specimen collection should always be considered in the interpretation of a laboratory result. Samples from early acute phase (i.e. those drawn before 3 days after rash onset) are more likely to test falsely IgM negative compared with those drawn 3 to 28 days after rash onset. For this reason, a second blood sample is indicated if the IgM serology results from an early acute phase sample are inconclusive or negative for measles, rubella, and parvovirus B19, and the person meets the clinical case definition for measles.

Sporadic cases, i.e. those cases with no epidemiologic link to a laboratory-confirmed case, nor having recent travel history to an area with known measles activity, must be laboratory-confirmed by measles virus isolation or have a demonstrated rise in IgG titre between acute and convalescent serum specimens. A false IgM positive result in a sporadic case is likely, even with highly specific IgM assays, due to the very low incidence of disease in Canada; therefore, IgM serology is not a reliable test in this particular situation.

### *Information on Virus Isolation and/or Culture*

Measles virus is present during the acute stage in throat and nasopharyngeal secretions for up to 4 days and is excreted in urine for at least 7 days after rash onset. Virus isolation should be attempted for all sporadic cases of measles and cases occurring early in a measles outbreak. With identification of a suspected measles case, the clinician should collect at least one appropriate specimen for virus isolation, in addition to blood for serology. If blood is not collected and the specimen obtained for virus isolation is required for laboratory confirmation of the case, the laboratory should be notified by phone. This should be documented on the requisition.

The isolation of measles virus may not always be successful. In an outbreak, specimens should be collected from several cases to increase the success of virus isolation and subsequent genotyping.

Direct detection of measles virus RNA in clinical specimens by reverse transcription-polymerase chain reaction is a potentially rapid option which can be used in addition to virus isolation. Contact the Viral Exanthemata Laboratory, LCDC, for procedural details.

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## The Public-Health Role

**The provincial or territorial epidemiologist should ensure that the following policy and procedures are in place locally.**

1. A method for informing and updating all physicians and local public-health authorities regarding the issues surrounding measles surveillance, specifically as follows:
  - Any case of rash illness meeting the clinical case definition for measles must be immediately reported to the local public-health authority.
  - Laboratory confirmation is essential for all cases of measles where there is no epidemiologic link to a confirmed case. A blood sample of at least 3 mL (e.g. full pediatric tube) should be taken for laboratory diagnosis.
  - The blood sample should be taken early in the illness, ideally between 3 and 7 days after onset of rash (at most 28 days after rash onset). If blood drawn before 3 days after rash onset is found to be negative for specific IgM antibodies against measles, rubella, and parvovirus B19 in a person who meets the clinical case definition for measles, a second specimen is indicated.
  - A blood specimen from suspected measles cases will also be screened for parvovirus B19 and rubella specific IgM. Although the differential diagnosis of rash is not limited to these viral illnesses, these are the most likely to be confused with measles and therefore the most important to rule out in Canada.
  - In addition to the blood sample, a urine specimen and/or nasopharyngeal swab or throat swab should be obtained for virus isolation. A nasopharyngeal or throat swab should be obtained within 4 days of onset of the rash and/or approximately 50 mL to 100 mL of sterile urine should be obtained within 7 days of onset of rash.
  - For successful measles virus isolation, the specimen must be collected, as described in Appendix D, and transported on ice (4° C) to the laboratory as soon as possible. Specimens must be processed by the laboratory within 48 hours after collection. LCDC can perform the sample processing and virus isolation should this be necessary.
  - If blood is not collected or cannot be obtained from a suspected measles case and a specimen obtained for viral isolation is available, the laboratory should be notified by phone that virus isolation will be required for laboratory confirmation of the case. This should be documented on the requisition.
  - It is essential for successful epidemiologic analysis of suspected measles cases, that a contact name and phone number, as well as the following information, be submitted to the laboratory with all specimens (see Appendix C):
    - name of institution sending specimen
    - patient identifier
    - patient name
    - date of birth
    - sex
    - city, county (municipality)
    - date of fever onset
    - date of rash onset

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In addition, the public-health authority will need to know the following:

- whether case meets clinical case definition (see Appendix A)
  - number of doses of measles vaccine received
  - date of last measles vaccination
  - date of collection
2. A means for transporting and processing of blood specimens for the serologic screen of rash.
  3. A means for transporting and processing of specimens taken for virus isolation.
  4. A policy and mechanism for the laboratory to notify the local public-health authority of all positive measles and rubella results within 24 hours of the results becoming available.
  5. A policy and mechanism for the local public-health authority to notify the provincial public-health authority of all positive measles and rubella results within one business day of the results becoming available.

### The Laboratory's Role

**The provincial or territorial laboratory must ensure the following regarding serologic testing.**

1. All serum specimens submitted from suspected measles cases are also screened for specific IgM antibodies against parvovirus B19 and rubella, constituting a "red rash screen," with a turnaround time of no more than 72 hours. Ideally these tests would all be done in a single laboratory.
2. Health professionals are aware of the appropriate timing of specimen collection for serologic testing and know that at least 3 mL of blood (e.g. one full pediatric tube) should be collected.
3. Only kits recommended by WGMEC are used to test for measles specific IgM (see Appendix E). The United States Centers for Disease Control and Prevention (CDC) measles IgM capture method, which is used at LCDC, is also a sensitive and specific assay which can be used for measles IgM serology. Serum specimens testing positive for measles IgM by a non-recommended assay are to be forwarded to LCDC for testing by the IgM capture assay or re-tested using a recommended commercial assay.
4. The laboratory that performs the serologic tests reports all results to the physician who ordered the test on the day they become available, and all positive results for measles and rubella to the local public-health authority within 24 hours of the results becoming available.

**The provincial or territorial laboratory must ensure the following regarding virus isolation and/or culture.**

1. Laboratories are prepared to advise clients as to the optimal time for specimen collection and provide handling instructions, specifically as follows:
  - A nasopharyngeal or throat swab must be collected within 4 days after onset of the rash and a sterile urine specimen

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(approximately 50 mL to 100 mL) within 7 days of rash onset, in order to increase the probability of isolating the virus.

- All specimens should be placed on ice and immediately transported to the laboratory for proper processing (see Appendix D).
  - Specimens should be processed within 48 hours, or frozen at  $-70^{\circ}\text{C}$  in viral transport medium (see Appendix D).
2. A procedure is in place for transferring specimens to a laboratory where virus isolation can be undertaken, and clients and the local public-health authorities are made aware of this procedure. For those laboratories that do not undertake virus isolation, the LCDC Viral Exanthemata Laboratory can provide sample processing and virus isolation services.
  3. Procedures are in place to send measles virus isolates from the testing laboratory to LCDC for genotypic analysis. LCDC should ensure that results of the genotypic analysis are forwarded to the isolating laboratory, within one business day of receipt of the results. The isolating laboratory should provide the results to the local public-health authority and to the provincial public-health authority within one business day of LCDC's notification.
  4. If the viral isolation was required for laboratory confirmation of a case, (i.e. if no serology was carried out), then
    - the result of the viral isolation attempt is communicated to the physician, who ordered the test, on the day it becomes available; and
    - the laboratory reports all positive results to the local public-health authority within 24 hours of the results becoming available.

**The provincial or territorial laboratory should ensure the following regarding proficiency testing.**

Since it is important to maintain a high level of proficiency in providing measles diagnostic services across the country, LCDC will offer proficiency testing programs. All laboratories providing measles diagnostic services should participate in these programs.

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## Appendix A

### NATIONAL SURVEILLANCE MEASLES CASE DEFINITIONS

#### Confirmed Case

A confirmed case is characterized by one of the following:

- laboratory confirmation of infection\* in the absence of recent (1 to 14 days) immunization with measles-containing vaccine:
  - isolation of measles virus from an appropriate clinical specimen
  - or, significant rise in measles specific antibody titre between acute and convalescent sera
  - or, positive serologic test for measles IgM antibody using a recommended assay
- **or**, clinical measles in a person who is epidemiologically linked to a laboratory-confirmed case

#### Clinical Case

A clinical case is characterized by all of the following:

- fever  $\geq 38.3^{\circ}\text{C}$
- cough, coryza, or conjunctivitis
- generalized maculopapular rash for at least 3 days

---

\* A laboratory-confirmed case does not have to meet the clinical illness description.

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
**Considerations**

- Laboratory confirmation of measles assumes the absence of recent measles vaccination.
- All suspected measles cases must test negative for parvovirus B19 and rubella virus.
- Interpretation of laboratory results must always consider the timing of specimen collection.
- The final decision for determination of a true measles case must consider clinical, laboratory, and epidemiologic data.

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**Appendix C**

**REQUISITION FOR VIRAL EXANTHEMATA**

Canadian Science Centre for Human and Animal Health  
Viral Exanthemata, LCDC  
1015 Arlington Street  
Winnipeg, Man. R3E 3R2  
Telephone: 204-789-6080

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### MEASLES VIRUS ISOLATION

#### Background

Measles virus isolation is a definitive means of establishing the diagnosis of measles. Measles culture has become a more reliable tool with the use of B95-8 cell line, an Epstein-Barr virus transformed lymphoblastoid cell line known to be quite sensitive in isolating measles virus from clinical specimens. As progress is made toward elimination of measles in the Americas, it will be critical to examine virus isolates from as many outbreaks and sporadic cases as possible for strain surveillance and to identify the source of the virus. From this standpoint, viral isolation is important. Automated DNA sequencing techniques are now available for rapid genetic characterization of viral isolates. This, together with the existence of a database of nucleic acid sequence information, now makes it possible to identify the source of wild-type viruses and rapidly differentiate between wild-type and vaccine strains.

Virus isolation can take several days to weeks. Therefore, **IgM serology should always be the first priority** in laboratory diagnosis of measles. Specimens for virus isolation should be taken **at the same time** that blood is obtained, since a delay in collection will reduce the chance of isolating the virus. Nasopharyngeal, or throat swab, or urine specimens should not be substituted for blood which is required for serologic diagnosis.

#### Protocols for Isolation of Measles Virus\*

Specimens for virus isolation should be obtained as soon as possible after the onset of rash. Always collect a urine specimen and, if possible, attempt to collect a nasopharyngeal or throat specimen.

#### Specimen Collection and Processing

##### *Respiratory Specimens (nasopharyngeal or throat swabs)*

1. Obtain the sample as soon as possible after onset of rash but no later than 4 days afterwards.
2. Use sterile swabs to obtain nasopharyngeal or throat specimens. The virus is cell-associated, so attempt to swab the throat and nasal passages to collect epithelial cells. Place the swabs in a tube containing 2 mL to 3 mL of VTM (viral transport medium: PBS or suitable isotonic solution such as Hank's BBS containing antibiotics [100 units/mL penicillin, 100 µg/mL streptomycin], and either 2% fetal bovine serum [FBS] or 0.5% gelatin).
3. Keep all specimens under refrigeration (4° C), and ship as soon as possible on wet ice to the laboratory for processing.

\* Original protocols, except for the shell vial method, are based on the method of Dr. William Bellini, Measles Section, CDC.

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4. In the laboratory, remove swabs from the VTM after allowing at least an hour for elution of the virus. Centrifuge the eluant at 2500 x g for 15 minutes at 4° C. Resuspend the pellet in 1 mL of tissue culture medium. If possible, save the supernatant in a separate tube. Either proceed directly with virus isolation (specimens may be held at 4° C for up to 24 hours), or freeze the samples at -70° C and ship on dry ice to the appropriate laboratory for virus isolation.

#### Urine Specimens

1. Collect 50 mL to 100 mL of urine within 7 days after rash onset.
2. Keep urine at 4° C and process within 48 hours at the latest.
3. Centrifuge 50 mL urine at 2500 x g for 15 minutes at 4° C to pellet the sediment. Resuspend the sediment in 2 mL VTM or any cell culture medium (e.g. DMEM, EMEM, RPMI plus antibiotics). Either proceed directly with virus isolation or freeze the samples at -70° C and ship on dry ice to the appropriate laboratory for virus isolation.

If centrifugation cannot be done, **do not freeze the urine sample**. The entire urine specimen should be held under refrigeration and shipped to the appropriate laboratory on wet ice. Seal the specimen container to prevent leakage.

#### Virus Isolation

##### B95-8 (B95-a) Cell Line

B95-8 is the preferred cell line for primary isolation of measles virus<sup>(1)</sup>, and this cell line is available from the American Type Culture Collection (No. CRL 1612). Note that **this cell line should be handled as an infectious cell line capable of yielding Epstein-Barr virus**.

The B95-8 cells grow lightly attached to the culture surface when grown in DMEM supplemented with 1 x antibiotics (100 units/mL penicillin, 100 mg/mL streptomycin), 0.25 µg/mL amphotericin (fungizone), and FBS. Cell growth is sustained by adding 5% to 10% FBS. Under these conditions, B95-8 cells grow as an adherent monolayer and are referred to as B95-a cells. FBS is used at a 2% concentration for cell maintenance during viral isolation. Grow cells in a moist CO<sub>2</sub> incubator at 37° C. Cell stocks can be stored frozen at -70° C using standard cryoprotection medium (50% FBS, 10% DMSO, 40% DMEM).

The B95-a cells can be passaged by briefly treating the cell monolayers with 0.05% trypsin-EDTA to release cells from the tissue culture surface. Be careful not to over-trypsinize. Neutralize trypsin by adding DMEM containing 10% FBS. Usually the cells from a single monolayer culture can be split in a ratio of 1:3. The cells tend to become "floaters" growing in clumps suspended in the medium as the cell density increases. These cells are viable and can be passaged by gentle pipetting to break up the clumps then replating to a lower cell density.

The B95-a cells can be transported at room temperature in a 75 cm<sup>2</sup> or 25 cm<sup>2</sup> tissue culture flask with additional medium added to help keep cells attached. After shipment, look at the cell sheet. If many cells are free-floating, a light spin of the medium will recover cells, which can be added back to the flask or to another flask for

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passage. For maintenance, add 30mL to 50 mL of DMEM with 2% FBS to a 75 cm<sup>2</sup> flask.

B95-a cells infected with measles virus can show syncytium formation and giant-cell cytopathic effect (CPE) as early as 48 hours after inoculation. However, isolation-attempt cultures should be followed for 7 to 8 days with two to three blind passages before isolation of measles virus is ruled out.

#### *Inoculation of Specimens for Measles Isolation (Flask Method – CDC)*

1. Passage cells, split into 25 cm<sup>2</sup> flasks at 1:3 or 1:4, and incubate 24 to 48 hours. Cells should be at 75% to 85% confluency when the specimen is inoculated.
2. For inoculation, decant medium, add 1.0 mL to 1.5 mL DMEM with 2 x antibiotics and add 0.1 mL to 1.0 mL of specimen, depending on the concentration. Save all leftover specimens as they may still be useful for PCR analysis.
3. Incubate at 37° C for 1 hour.
4. Decant medium into bleach water; replace medium with 5mL to 10 mL of DMEM containing 2% FBS and 2 x antibiotics.
5. Change medium every 3 to 4 days, and passage cells by splitting at 1:3 every 7 to 9 days. Check for CPE daily.
6. Attempt at least three blind passages before terminating isolation attempt.
7. If CPE is visible, continue to feed the cells until the CPE becomes extensive. It may be necessary to passage the cells one more time to allow the CPE to progress. When CPE is maximal, pellet cells, resuspend the pellet in 1 mL of DMEM, and freeze at –70° C.
8. If a culture becomes contaminated, the original specimen can be diluted in DMEM and passed through a 0.45 µm nitrocellulose filter and used to inoculate fresh B95-a cells.
9. Confirmation of culture can be achieved by using fluorescent antibody staining or PCR.

**Please remember to save some of the original clinical specimen. This material can be used for a second isolation attempt if problems occur with the first, and it can provide a specimen for direct RT-PCR analysis. (Contact the Viral Exanthemata Laboratory, LCDC, for details.)**

#### *Inoculation of Specimens for Virus Isolation (Shell Vial Method)*

The shell vial method has been found to be quite convenient and highly successful by some laboratories in Canada. This method involves the following:

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### Maintenance of B95-8 cells

1. B95-8 cells are routinely maintained in a 25 cm<sup>2</sup> flask in an upright position. The cells, contained in 10 mL RPMI growth medium, settle to the bottom of the flask.
2. Without disturbing the cells, each week carefully remove 9 mL of growth medium and discard, leaving 1 mL of cell suspension. Add 9 mL of fresh growth medium (formula below). If additional cells are required, each flask can be split 1:3.
3. Trypsinization is not required because these cells do not attach. Reduction in the density of the cells can be accomplished by harvesting the cells through centrifugation and resuspending them in the fresh medium to the required density.
4. These cells can be maintained at -70° C in an appropriate cryoprotectant and easily retrieved when required.

### Growth medium for B95-8 cells

- 500 mL RPMI 1640 (Gibco Cat. No. 21870-076 without glutamine).
- Aseptically remove and save 56 mL of the medium for other use.
- Add 50 mL FBS (heat inactivated) to the remaining 444 mL medium.
- Add 1 mL penicillin-gentamicin solution (final concentration: penicillin 100 units/mL, gentamicin 10 µg/mL).
- Add 5 mL 300 mM glutamine.

### Inoculation

1. Inoculate two shell vials (containing coverslips) per specimen with 10<sup>6</sup> cells in 0.2 mL growth medium and 0.2 mL of specimen.
2. Centrifuge 1,000 x g for 45 minutes at room temperature.
3. Following centrifugation, remove inoculum and add 1 mL maintenance medium containing 2% FBS.
4. Incubate at 37° C and observe for syncytium formation at 48 hours and daily for up to 6 days.
5. When CPE is observed (at least 50% cells), or after 6 days incubation, centrifuge one vial at 900 x g for 10 minutes.

#### **If viral isolation is needed to confirm the case, follow steps 6 and 7.**

6. Aspirate the maintenance medium from the vial, fix by adding 1 mL acetone and allow to sit for 15 minutes. Remove the acetone and wash with PBS to remove residual acetone.
7. Aspirate PBS and stain the coverslip (within the vial) with measles immunofluorescent stain according to the manufacturer's instructions (Light Diagnostics – Chemicon; distributed in Canada by Bio/Can).

#### **Otherwise inoculate fresh B95-8 cells in 75 cm<sup>2</sup> flask with virus until CPE is observed, and follow steps 7 and 8 of the Flask Method protocol.**

8. If CPE is not observed after 6 days, the monolayer in the second shell vial can be scraped and passed to a fresh vial(s) for further incubation.

**Remember to save some of the original clinical specimen. This material can be used for a second isolation attempt if**

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problems occur with the first, as well as provide a specimen for direct RT-PCR analysis. (Contact the Viral Exanthemata Laboratory, LCDC, for details.)

#### Shipping of Measles Virus Isolates

Infected cells can be pelleted, resuspended in a small volume of DMEM and frozen at  $-70^{\circ}\text{C}$  before shipping on dry ice.

The appropriate Transport Canada regulations for the safe shipment of infectious material and dangerous goods must always be followed.

Please contact LCDC concerning shipment of sera, clinical specimens, or virus isolates.

**Viral Exanthemata Laboratory**  
**Canadian Science Centre for Human and Animal Health**  
**Laboratory Centre for Disease Control**  
**1015 Arlington Street**  
**Winnipeg, Man. R3E 3R2**

**Head: Dr. Graham Tipples**  
**Telephone: 204-789-6080**  
**Fax: 204-789-5009**  
**E-mail: graham\_tipples@hc-sc.gc.ca**

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## Appendix E

### CURRENT RECOMMENDATIONS REGARDING COMMERCIAL MEASLES IgM SEROLOGY TEST KITS

The laboratory subcommittee of the WGMEC recently completed a study on the relative performance of indirect measles IgM enzyme immunoassays (EIAs) and IgM capture assays using a total of 308 serum samples from confirmed measles outbreak cases and 454 samples from non-measles cases<sup>(2)</sup>. This evaluation included the Behring, Clark, Gull, and PanBio EIAs and a commercial IgM capture assay (Light Diagnostics) along with the CDC measles IgM capture assay. Overall, the Gull and Behring EIAs (sensitivity: 88.3%, 88.6%; specificity: 99.6%, 96.7%; positive predictive value: 99.6%, 97.8%; negative predictive value: 93.0%, 94.6%, respectively) were found to be the best commercial kits for measles IgM serology. The CDC measles IgM capture assay was also found to be highly sensitive and specific. The CDC measles IgM capture procedure is available at the LCDC. Based on this data, the WGMEC currently recommends the use of either the Gull or Behring commercial EIAs for measles IgM serology.

The timing of the specimen collection was an important criterion for IgM seropositivity regardless of the test used, and the best seropositivity rates in the range of 92% to 100% were observed with samples collected 6 to 14 days after the onset of symptoms. Sera positive for other IgM markers such as rubella, parvovirus B19, EBV, and HHV-6, etc. gave false positive reactions in some of the assays including the capture assays.

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	<b>Behring</b>	<b>Clark</b>	<b>Gull</b>	<b>PanBio</b>	<b>Light Diagnostics</b>	<b>CDC</b>
Sensitivity <sup>(a)</sup> Sensibilité <sup>(a)</sup>	88.6% (85.1-92.1) <sup>(b)</sup>	82.8% (76.6-87.0)	88.3% (84.7-91.9)	83.1% (78.9-87.3)	92.2% (89.2-95.2)	87.0% (82.3-91.7)
Specificity <sup>(c)</sup> Spécificité <sup>(c)</sup>	96.7% (95.1-98.3)	97.1% (95.6-98.6)	99.6% (99.0-100)	86.6% (83.5-89.7)	86.6% (83.5-89.7)	94.8% (92.4-96.8)
Positive predictive value Valeur prévisionnelle positive	97.8% (96.1-99.5)	95.9% (93.5-98.3)	99.6% (98.9-100)	97.7% (95.9-99.5)	88.2% (84.7-91.7)	96.1% (93.3-98.9)
Negative predictive value Valeur prévisionnelle négative	94.6% (92.5-96.7)	90.2% (87.6-92.8)	93.0% (90.7-95.3)	91.4% (88.8-94.0)	95.9% (94.0-97.8)	95.7% (93.5-97.5)

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## References

1. Kobune F, Sakata H, Sugiura A. *Marmoset lymphoblastoid cells as a sensitive host for isolation of measles virus*. J Virol 1990;64:700-05.
2. Ratnam S, Tipples G, Head C et al. *Performance of indirect immunoglobulin M (IgM) serology tests and IgM capture assays for laboratory diagnosis of measles*. J Clin Microbiol. In press.

## Références

1. Kobune F, Sakata H, Sugiura A. *Marmoset lymphoblastoid cells as a sensitive host for isolation of measles virus*. J Virol 1990;64:700-05.
2. Ratnam S, Tipples G, Head C et coll. *Performance of indirect immunoglobulin M (IgM) serology tests and IgM capture assays for laboratory diagnosis of measles*. J Clin Microbiol (Sous presse).

***Our mission is to help the people of Canada maintain and improve their health.***

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*Santé Canada*

The Canada Communicable Disease Report (CCDR) presents current information on infectious and other diseases for surveillance purposes and is available through subscription. Many of the articles contain preliminary information and further confirmation may be obtained from the sources quoted. Health Canada does not assume responsibility for accuracy or authenticity. Contributions are welcome (in the official language of your choice) from anyone working in the health field and will not preclude publication elsewhere.

Scientific Advisors	Dr. John Spika	(613) 957-4243
	Dr. Fraser Ashton	(613) 957-1329
Editor-in-Chief	Eleanor Paulson	(613) 957-1788
Assistant Editor	Nicole Beaudoin	(613) 957-0841
Desktop Publishing	Francine Boucher	

Submissions to the CCDR should be sent to the Editor-in-Chief, Laboratory Centre for Disease Control, Tunney's Pasture, Address Locator 0602C2, Ottawa, Ontario K1A 0L2.

To subscribe to this publication, please contact:

Canadian Medical Association	Tel. No.:	(613) 731-8610 Ext. 2307
Member Service Centre		or (888) 855-2555
1867 Alta Vista Drive	FAX:	(613) 236-8864
Ottawa, ON Canada K1G 3Y6		

Annual subscription: \$83.00 (plus applicable taxes) in Canada; \$109 (U.S.) outside Canada.

© Minister of Health 1999 (On-line) ISSN 1481-8531  
Publications Mail Agreement No. 1437887

This publication can also be accessed electronically via Internet using a Web browser at <<http://www.hc-sc.gc.ca/hpb/lcdc>>. It can also be accessed at any time from any fax machine using LCDC's FAXlink Service by calling 1-613-941-3900.

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