

**Positive Result Thresholds**  
for the Smith Detection (Barringer) IONSCAN

Subject to the specific equipment being used, analysis results at or above the levels identified below will be considered to represent a positive (+) result.

**IONSCAN (Barringer) System**

<b><u>Substance</u></b>	<b><u>Alarm threshold</u></b>
Cocaine (includes crack)	<b>500</b>
Heroin	<b>100</b>
Amphetamine (uppers)	<b>500</b>
Methamphetamine (uppers)	<b>500</b>
MDA (ecstasy)	<b>300</b>
MDMA (ecstasy)	<b>300</b>
MDEA (ecstasy)	<b>300</b>
THC (active component of hash/marijuana)	<b>50</b>
LSD (acid)	<b>100</b>
PCP	<b>500</b>
Hash/marijuana (grass, pot)	<b>Any alarm</b>
Opium	<b>Any alarm</b>

**Résultat positif – Seuils minimaux** applicables  
à l'IONSCAN de Smith Detection (Barringer)

Selon le type d'appareil utilisé, tout résultat d'analyse qui correspond ou est supérieur aux seuils indiqués ci-dessous sera considéré comme un résultat positif (+).

**Système IONSCAN (Barringer)**

<b><u>Substance</u></b>	<b><u>Seuil d'alarme</u></b>
Cocaïne (y compris le crack)	<b>500</b>
Héroïne	<b>100</b>
Amphétamine (speed)	<b>500</b>
Méthamphétamine (speed)	<b>500</b>
MDA (ecstasy)	<b>300</b>
MDMA (ecstasy)	<b>300</b>
MDEA (ecstasy)	<b>300</b>
THC (composante active du hasch/marijuana)	<b>50</b>
LSD (acide)	<b>100</b>
PCP	<b>500</b>
Hasch/marijuana (herbe, pot)	<b>N'importe quelle alarme</b>
Opium	<b>N'importe quelle alarme</b>